

Baso

CYTOLOGY STAIN SERIES

PAPANICOLAOU STAIN

Papanicolaou Staining is used to differentiate cells in smear preparations in various bodily secretions; the specimens can be gynaecological smears (Pap Smears), liquid based cytology preparations, sputum, brushings, washings, urine, cerebrospinal, abdominal, pleural, synovial and seminal fluids, fine needle aspiration material, tumor touch samples, or other materials containing cells.

Pap Staining involves five dyes in three solutions, a nuclear stain and two counter stains:

Nuclear Stain: Harris Hematoxylin which is used to stain cell nuclei.

Counter Stain OG6: Orange G dye is used to stain keratin. Its original role was to stain the small cells of keratinizing squamous cells present in sputum.

Counter Stain EA50: EA50 (Eosin Azure) counter stain, comprising three (3) dyes; Eosin Y, Light Green SF Yellowish and Bismark Brown.

- Eosin Y stains the superficial epithelial squamous cells, nucleoli, cilia, and red blood cells
- Light Green yellowish stains the cytoplasm of other cells. Including non-keratinized squamous cells.
- Bismark Brown Y which is used to balance colour intensity

Hazards Identified:

Contains:

Methanol <5%

EC No.: 200-659-6



FLAMMABLE



SYSTEMATIC
HEALTH HAZARDS



TOXIC



Available in:

Papanicolaou Stain Kit (Hematoxylin, EA-50, OG6) - 3 x 250ml/kit

Papanicolaou Stain (EA-50) - 1 Liter

Papanicolaou Stain (OG-6) - 1 Liter



TruLaboratories
Corporation

+632 82-TRULAB (878522)

+63917.878.2288 • +63925.878.2288

www.trulaboratories.com